

COVER GOWN TECHNICAL SPECIFICATION

Material	100% PPSB with PE coating
Style	Disposable Coverall
Application	Medical & Health
Characteristic	<ul style="list-style-type: none">- Breathable and impervious- Liquid Resistant- Hood design provides overall protection- Suitable for health care workers and visitors, helps to prevent contamination with particulates.- Full coverage, breathable and comfortable- Waist and neck ties, Elastic or Knit Cuffs
Size	L-XL EXPANDIBLE
Weight	40 gsm
Certification	ISO4920

