

HAIR CAP TECHNICAL SPECIFICATION

Material	100% PPSB with PE coating
Style	Disposable Coverall
Application	Medical & Health
Characteristic	- Breathable and impervious - Liquid Resistant - Suitable for health care workers and visitors, helps to prevent contamination with particulates. - Comfortable
Weight	40 gsm
Color	White
Certification	ISO4920

Ensure a clean workspace with hair cap. With its flexible, pleated design, this hair cap easily stretches over the worker's head to keep hair away from the face and prevent loose strands from contaminating sensitive materials and products. It comes in a white color for a clean appearance, and they're disposable for fresh, hygienic use every time. Hair cap is made of durable polypropylene and ideal for use in food service and health facilities.

- Flexible, pleated design that stretches over the worker's head
- White color for a clean appearance
- Disposable design for fresh, hygienic use every time
- Polypropylene material for use in food service and health facilities

