

PROTECTIVE COVERALL (LAMINATED) TECHNICAL SPECIFICATION

Material	100% PPSB with PE coating
Style	Disposable Coverall
Application	Medical & Health
Characteristic	<ul style="list-style-type: none"> - Breathable and impervious - Synthetic blood penetration resistant test passed - Hood design provides overall protection - With elastic wrist and open ankle
Size	L-XL EXPANDIBLE
Weight	55 gsm
Color	White
Certification	ISO4920



Anatomical design and fit for optimal wear, coverage and comfort

Two-way zipper for on/off convenience



Storm flap on all coverall to help provide additional protection

Increased comfort and mobility